



DECEMBER 11-12, 2025
۲۱ و ۲۰ آذر ماه ۱۴۰۴



دومین کنگره
پلاسما پزشکی ایران

The 2nd Congress on Plasma Medicine

دبیرخانه دائمی کنگره
پلاسما پزشکی ایران
www.plasmamedsym.ir



Unprecedented Eradication of Biofilms on Flexible Endoscopes Using Cold Atmospheric Plasma-Aerosol

Mohammad Mahdi Ajalli

Department of Medicine, Zanjan University of Medical Sciences, Zanjan, Iran

OPEN ACCESS

Citation:

Ajalli MM. Unprecedented Eradication of Biofilms on Flexible Endoscopes Using Cold Atmospheric Plasma-Aerosol. *Iran Biomed J. Supplementary* (2026): 50.



This article is licensed under a Creative Commons Attribution-NonDerivatives 4.0 International License.

Keywords: Biofilm inactivation, Cold atmospheric plasma, Endoscope sterilization, Low-temperature disinfection, Reactive species

ABSTRACT

Introduction: Sterilization of medical equipment is critical to prevent healthcare-associated infections, yet conventional methods often fail against resilient biofilms on heat-sensitive devices. However, effective low-temperature alternatives that eradicate biofilms without damaging materials remain elusive. In this study, we demonstrate the unprecedented efficacy of cold atmospheric plasma-aerosol (CAP-A) in sterilizing biofilm-contaminated endoscopes.

Materials and Methods: We exposed *Pseudomonas aeruginosa* and *Staphylococcus aureus* biofilms on flexible endoscopes to CAP-A generated from a helium-oxygen mixture at room temperature. Treatment durations ranged from 5 to 15 minutes, followed by viability assessments via colony-forming unit counts and confocal microscopy.

Results and Discussion: CAP-A treatment reduced biofilm viability by over 6 log units (99.9999% inactivation) within 10 minutes, compared to less than 2 log units for hydrogen peroxide vapor ($n = 12$; $p < 0.001$). Reactive oxygen and nitrogen species penetrated biofilms up to 50 μm deep, resolving structures at 1 μm resolution via microscopy. No material degradation occurred, with surface integrity preserved (>99% similarity to controls). Against multi-species biofilms, inactivation reached 7 log units, establishing superior performance over prevailing low-temperature methods.

Conclusion: The approach suggested in our study establishes CAP-A as a transformative sterilization technique for thermolabile medical devices, paving the way for reduced infection rates in clinical settings.

Corresponding Author: Mohammad Mahdi Ajalli

Department of Medicine, Zanjan University of Medical Sciences, Zanjan, Iran; E-mail: Zums.mma@gmail.com



Iranian Biomedical Journal Supplementary (February 2026): 50